			AB2-5 AB4-8	AB2-5 (gray lens) AB2-7	AC2-5	4C-A	SP4-10 SP10-16	SP6-12	PA2-5P	PAA-7 PAG-8 PA2-5	IC5-9	RAB4-8P RAB2-5	RSP5-12 RSP6-12 RSP6-16	RIC5-9 RIC5-9H RIC5-9W	RRE6-10	SCW2.0	PCW4.0	RAB4-8L RAB2-5L	RNA5-9		
GE Ultrasound Contact Gel	:Gel	:GE Accessories	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-
Aquasonics 100	:Gel	:Parker Laboratory Inc.	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-
HAIYIN Gel	:Gel	:Wuxi Huaili Industrial Company	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-
PROJELLY	:Gel	:JEX Inc.	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-
Sonogel	:Gel	:Sonogel Vertriebs GmbH	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-
Scan	:Gel	:Parker Laboratory Inc.	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-
Shanghai A	:Gel	:Shanghai She Li Yiliao Qxie Co.Ltd	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-
Shanghai B	:Gel	:Shanghai She Li Yiliao Qxie Co.Ltd	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-
Wavelength	:Gel	:National Therapy Products Inc.	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	-
Incidin Foam	:20%IPA, 10%Ethyl Alcohol	:Ecolab	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 5 minutes
T-Spray	:Quaternary Ammoniums	:Pharmaceutical Innovations	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 20 minutes
T-Spray II	:Quaternary Ammoniums	:Pharmaceutical Innovations	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 20 minutes
Transeptic Spray	:70% IPA	:Parker Laboratory Inc.	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
Asepti-Wipes	:8%IPA, <1%Quat	:Ecolab Inc.	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
Asepti-Wipes II	:14.85%IPA	:Ecolab Inc.	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
Cleanisep wipes	:Quaternary Ammoniums	:Dr.Schumacher GmbH	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
Septiwipes	:Quaternary Ammoniums	:Dr.Schumacher GmbH	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
CaviWipes	:14.3%IPA	:Metrex	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
MetriWipes	:8%IPA, <1%Quat	:Metrex	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
NUCLEAN II	:Isopropanol / Ethanol	:Alkapharm	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
Sani-Cloth Active	:Quaternary Ammoniums	:PDI	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
Sani-Cloth Germicidal	:8%IPA, <1%Quat	:PDI	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
Sani-Cloth HB	:Quaternary Ammoniums	:PDI	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 30 minutes
PeraSafe	:Peracetic Acid	:Antec International	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
Alkazyme	:Proteolytic Enzymes	:Alkapharm	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 10 minutes
Cidex	:Glutaraldehyde	:Johnson & Johnson	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 50 minutes
Cidex OPA	:Ortho-phthlaldehyde	:Johnson & Johnson	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
Cidex Plus	:Glutaraldehyde	:Johnson & Johnson	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 25 minutes
Gigasept	:Succindialdehyde, formaldehyde	:Schülke & Mayr	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 30 minutes
Gigasept AF	:Quaternary Ammoniums	:Schülke & Mayr	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 30 minutes
Gigasept FF	:Succindialdehyde	:Schülke & Mayr	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 30 minutes
Klenzyme	:Proteolytic Enzymes	:STERIS	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 10 minutes
Metricide 14	:Glutaraldehyde	:Metrex	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 50 minutes
Metricide 28	:Glutaraldehyde	:Metrex	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 50 minutes
Metricide 30	:Glutaraldehyde	:Metrex	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 25 minutes
MetriCide OPA Plus	:Ortho-phthlaldehyde	:Metrex	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 15 minutes
Sekusept Plus	:Nonionic surfactant	:Henkel-Ecolab	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 30 minutes
Ster-Bac Blu	:Quaternary Ammoniums	:Ecolab Inc.	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 50 minutes
Wavicide-01	:Glutaraldehyde	:MCC	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	OK	< 50 minutes

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