			RSP6-16-RS	RIC5-9-RS RIC5-9W-RS	RAB2-5-RS	RAB4-8-RS	RNA5-9-RS	AB2-7-RS	SP102-16-RS	
	GE Ultrasound Contact Gel	Gel	GE Accessories	OK	OK	OK	OK	OK	OK	—
	Aquasonics 100	Gel	Parker Laboratory Inc.	OK	OK	OK	OK	OK	OK	—
	Haiyin	Gel	Wuxi Huaili Industrial Company	⊗	⊗	⊗	⊗	⊗	OK	—
	PROJELLY	Gel	JEX Inc.	⊗	⊗	⊗	⊗	⊗	OK	—
	Sonogel	Gel	Sonogel Vertriebs GmbH	OK	OK	OK	OK	OK	OK	—
	Scan	Gel	Parker Laboratory Inc.	OK	OK	OK	OK	OK	OK	—
	Incidin Foam	20%IPA, 10%Ethyl Alcohol	Ecolab	⊗	OK	⊗	⊗	⊗	⊗	< 5 minutes
	T-Spray	Quaternary Ammoniums	Pharmaceutical Innovations	OK	OK	OK	OK	OK	⊗	< 20 minutes
	T-Spray II	Quaternary Ammoniums	Pharmaceutical Innovations	OK	OK	OK	OK	OK	⊗	< 20 minutes
	Transeptic Spray	70% IPA	Parker Laboratory Inc.	⊗	OK	OK	OK	⊗	OK	< 15 minutes
	Asepti-Wipes	8%IPA, <1%Quat	Ecolab Inc.	⊗	OK	OK	OK	OK	⊗	< 15 minutes
	Septiwipes	Quaternary Ammoniums	Dr.Schumacher GmbH	OK	OK	OK	OK	⊗	⊗	< 15 minutes
	Cleanisept	Quaternary Ammoniums	Dr.Schumacher GmbH	OK	OK	OK	OK	⊗	⊗	< 15 minutes
	MetriWipes	8%IPA, <1%Quat	Metrex	⊗	OK	OK	OK	⊗	⊗	< 15 minutes
	NUCLEAN II	Isopropanol / Ethanol	Alkapharm	⊗	OK	⊗	⊗	⊗	⊗	< 15 minutes
	Sani-Cloth Active	Quaternary Ammoniums	PDI	OK	⊗	⊗	⊗	⊗	⊗	< 15 minutes
	Sani-Cloth Germicidal	8%IPA, <1%Quat	PDI	⊗	OK	OK	OK	OK	⊗	< 15 minutes
	Sani-Cloth HB	Quaternary Ammoniums	PDI	OK	OK	OK	OK	OK	OK	< 30 minutes
	PeraSafe	Peracetic Acid	Antec International	OK	OK	OK	OK	OK	⊗	< 15 minutes
	Alkazyme	Proteolytic Enzymes	Alkapharm	⊗	⊗	⊗	⊗	⊗	⊗	< 10 minutes
	Cidex	Glutaraldehyde	Johnson & Johnson	OK	OK	OK	OK	OK	OK	< 50 minutes
	Cidex OPA	Ortho-phthalaldehyde	Johnson & Johnson	OK	OK	OK	OK	OK	OK	< 15 minutes
	Cidex Plus	Glutaraldehyde	Johnson & Johnson	⊗	OK	OK	⊗	OK	⊗	< 25 minutes
	Gigasept	Succindialdehyde, formaldehyde	Schülke & Mayr	OK	OK	⊗	⊗	OK	⊗	< 30 minutes
	Gigasept AF	Quaternary Ammoniums	Schülke & Mayr	⊗	OK	OK	OK	⊗	⊗	< 30minutes
	Gigasept FF	Succindialdehyde	Schülke & Mayr	OK	OK	OK	OK	⊗	⊗	< 30minutes
	Klenzyme	Proteolytic Enzymes	STERIS	⊗	OK	OK	OK	⊗	⊗	< 10 minutes
	Metricide 14	Glutaraldehyde	Metrex	OK	OK	OK	OK	OK	OK	< 50 minutes
	Metricide 28	Glutaraldehyde	Metrex	⊗	⊗	⊗	⊗	⊗	⊗	< 50 minutes
	Metricide 30	Glutaraldehyde	Metrex	⊗	⊗	⊗	⊗	⊗	⊗	< 25 minutes
	Metricide OPA Plus	Ortho-phthalaldehyde	Metrex	OK	OK	OK	OK	OK	OK	< 15 minutes
	Sekusept Plus	Nonionic surfactant	Henkel-Ecolab	⊗	OK	⊗	⊗	⊗	⊗	< 30 minutes
	Ster-Bac Blu	Quaternary Ammoniums	Ecolab Inc.	OK	OK	OK	OK	OK	OK	< 50 minutes
	Wavicide-01	Glutaraldehyde	MCC	OK	OK	OK	OK	OK	OK	< 50 minutes

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