

# GE CRITIKON Blood Pressure Cuff Adapters and Connectors

Total solutions for non-invasive blood pressure.  
Connecting intelligence and care.



# GE's connection system is simple and easy to implement.

GE Healthcare provides cost-effective standardization and inventory management solutions. By standardizing blood pressure (BP) connection systems, you create a universal connection to allow the same type of BP cuff connection to work on a wide variety of non-invasive blood pressure devices. By consolidating all of your BP cuff purchases to one vendor, your Purchasing Department will have fewer purchase orders and invoices. And a single point of contact for BP cuffs is created, providing you with a dedicated implementation team.

Products included in this catalog are not made with natural rubber latex.

## BP Cuff Adapters

### Connecting BP Cuff to 2 Male Slip Luers



**Y, Unfitted to 2 Male Slip Luers**  
Compatibility: Unfitted  
**Part # 330062 (Qty: 10)**



**Y, 1 Male Bayonet to 2 Male Slip Luers**  
Compatibility: Philips/HP,<sup>®</sup> Siemens,<sup>®</sup> Datascope,<sup>®</sup> Colin<sup>®</sup>  
**Part # 330071 (Qty: 10)**



**Y, 1 Female Locking Luer to 2 Male Slip Luers**  
Compatibility: Spacelabs,<sup>®</sup> Datascope, Colin  
**Part # 330065 (Qty: 10)**



**Y, 1 Female Submin to 2 Male Slip Luers**  
Compatibility: Protocol<sup>®</sup>  
**Part # 330066 (Qty: 10)**



**X, 2 Female Submins to 2 Male Slip Luers**  
Compatibility: GE, GE Marquette\*  
**Part # 300939 (Qty: 10)**



**Y, 1 Male Submin to 2 Male Slip Luers**  
Compatibility: Nihon Kohden<sup>®</sup>  
**Part # 330069 (Qty: 10)**



**X, 2 Mated Pair Submins to 2 Male Slip Luers**  
Compatibility: GE Datex-Ohmeda\*  
**Part # 330063 (Qty: 10)**



**X, 2 Male Bayonets to 2 Male Slip Luers**  
Compatibility: User Modified  
**Part # 330067 (Qty: 10)**

## Connecting BP Cuff to Female Slip Luer



**Y, 1 Male Slip Luer to  
2 Female Slip Luers-Neonatal**  
Compatibility: Philips/HP, Siemens,  
Datascope, Spacelabs  
**Part # 330049 (Qty: 1)**



**1 Male Micro Plug to  
1 Female Locking Luer-Neonatal**  
Compatibility: Spacelabs  
**Part # 330088 (Qty: 10)**



**Y, 1 Male Micro Plug to 2 Female  
Locking Luers-Neonatal**  
Compatibility: Spacelabs  
**Part # 330078 (Qty: 10)**

## Connecting BP Cuff to 2 Female Screw Connectors



**Y, 1 Unfitted to  
2 Female Screw Connectors**  
Compatibility: Unfitted  
**Part # 330048 (Qty: 1)**



**Y, 1 Male Bayonet to  
2 Female Screw Connectors**  
Compatibility: Philips/HP, Siemens,  
Datascope, Colin  
**Part # 330046 (Qty: 1)**

## Connecting BP Cuff to 2 Unfitted



**Y, 1 Unfitted to 2 Unfitted**  
Compatibility: Unfitted  
**Part # 300829 (Qty: 1)**

# BP Cuff Adapters

## Connecting BP Cuff to 2 Male Submins



**Y, 1 Unfitted to 2 Male Submins**  
Compatibility: Unfitted  
**Part # 330084 (Qty: 10)**



**Y, 1 Male Bayonet to 2 Male Submins**  
Compatibility: Philips/HP, Siemens, Datascope, Spacelabs  
**Part # 330080 (Qty: 10)**



**Y, 1 Female Locking Luer to 2 Male Submins**  
Compatibility: Legacy Spacelabs, Datascope, Colin  
**Part # 330081 (Qty: 10)**



**Y, 1 Female Submin to 2 Male Submins**  
Compatibility: Protocol  
**Part # 330093 (Qty: 10)**



**Y, 1 Male Submin to 2 Male Submins**  
Compatibility: Nihon Kohden  
**Part # 330083 (Qty: 10)**



**X, 2 Male Screw Connectors to 2 Male Submins**  
Compatibility: GE DINAMAP\*  
**Part # 330087 (Qty: 10)**



**X, 2 Male Bayonets to 2 Male Submins**  
Compatibility: MDE, Spacelabs  
**Part # 330086 (Qty: 10)**



**X, Mated Pair Submins to 2 Male Submins**  
Compatibility: GE Datex-Ohmeda  
**Part # 330085 (Qty: 10)**



# BP Cuff Adapters

## Connecting BP Cuff to NEW DINACLICK and Neo-Snap Connectors



**Y, 1 Male Bayonet to DINACLICK\* Female Recept**

Compatibility: Philips/HP, Siemens, Datascope, Colin

Part # 2059291-001 (Qty: 10)



**Y, 1 Female Submin to DINACLICK Female Recept**

Compatibility: Protocol

Part # 2059291-002 (Qty: 10)



**Y, T Male Submin to DINACLICK Female Recept**

Compatibility: Nihon Kohden

Part # 2059291-003 (Qty: 10)



**Y, 1 Unfitted to DINACLICK Female Recept**

Compatibility: Unfitted

Part # 2059291-004 (Qty: 10)



**X, 2 Male Bayonet to DINACLICK Female Recept**

Compatibility: User Modified

Part # 2059292-001 (Qty: 10)



**X, 2 Female Submin to DINACLICK Female Recept**

Compatibility: GE Marquette

Part # 2059292-002 (Qty: 10)



**X, 2 Mated Submin to DINACLICK Female Recept**

Compatibility: GE Datex-Ohmeda

Part # 2059292-003 (Qty: 10)



**X, 2 Male Screw to DINACLICK Female Recept**

Compatibility: DINAMAP

Part # 2059292-004 (Qty: 10)



**Y, 1 Male Slip Luer to Neo-Snap\* Male Plug**

Compatibility: Philips/HP, Siemens, Datascope, Spacelabs

Part # 2059290-001 (Qty: 10)



**Male Micro Plug to 2 Neo-Snap Male Plug**

Compatibility: Spacelabs

Part # 2059290-002 (Qty: 10)



**Male Micro Plug to 1 Neo-Snap Male Plug**

Compatibility: Spacelabs

Part # 2059297-001 (Qty: 10)

# BP Cuff Connectors



**5/32" ID Barbed to 1/8" ID Barbed, Plastic**  
 Compatibility: Airhose Connector  
**Part # 300671 (Qty: 10)**



**5/32" ID Barbed to 5/32" ID Barbed, Plastic**  
 Compatibility: Airhose Connector  
**Part # 300670 (Qty: 10)**



**Female Submin to 5/32" ID Barbed, Plastic**  
 Compatibility: GE, GE Marquette, GE Datex-Ohmeda  
**Part # 330092 (Qty: 10)**



**5/32" ID Barbed to Male Slip Luer, Plastic**  
 Compatibility: Female Locking Luer and Screw Connectors  
**Part # 300666 (Qty: 10)**



**Female Locking Luer to 5/32" ID Barbed, Plastic**  
 Compatibility: Spacelabs, Datascope, Colin  
**Part # 300668 (Qty: 10)**



**Male Screw to Male Slip Luer, Plastic**  
 Compatibility: Female Locking Luer and Screw Connectors  
**Part # 300662 (Qty: 10)**



**Male Screw to 5/32" ID Barbed, Plastic**  
 Compatibility: GE DINAMAP  
**Part # 300664 (Qty: 10)**



**5/32" ID Barbed to Female Bayonet, Metal**  
 Compatibility: Philips/HP, Siemens, Datascope, Colin  
**Part # 330060 (Qty: 10)**



**Male Bayonet to 3/16" ID Barbed, Metal**  
 Compatibility: Philips/HP, Siemens, Datascope, Colin  
**Part # 330059 (Qty: 10)**



**2T DINACLICK Female Recept to 1/8" ID tubing, dual lumen hoses**  
 Compatibility: DINACLICK Male Plug  
**Part # 2059294-001 (Qty: 10)**



**2-1T DINACLICK Female Recept for single lumen hoses**  
 Compatibility: DINACLICK Male Plug  
**Part # 2055630-001 1/8" tube (Qty: 10)**  
**Part # 2055630-002 5/32" tube (Qty: 10)**



**5/32" ID Barbed to Female Slip Luer**  
 Compatibility: Universal Neonatal  
**Part # 330068 (Qty: 10)**



**5/32" ID Barbed to Female Legacy Marquette Connector, Plastic**  
 Compatibility: GE Marquette  
**Part # 330057 (Qty: 10)**



**2T Neo-Snap Male Plug**  
 Compatibility: Neo-Snap Female Care hose  
**Part # 2059293-001 (Qty: 10)**



**Male Submin to  
1/8" ID Barbed, Plastic**  
Compatibility: GE, GE Marquette,  
GE Datex-Ohmeda  
**Part # 330091 (Qty: 10)**



**Male Submin to  
5/32" ID Barbed, Plastic**  
Compatibility: GE, GE Marquette,  
GE Datex-Ohmeda  
**Part # 330090 (Qty: 10)**



**Mated Pair Submin to  
1/8" ID Barbed, Plastic**  
Compatibility: GE Datex-Ohmeda  
**Part # 330058 (Qty: 10)**



**Female Submin to  
1/8" ID Barbed, Plastic**  
Compatibility: GE,  
GE Marquette, Protocol  
**Part # 330064 (Qty: 10)**



**1/8" ID Barbed to  
Female Screw, Metal**  
Compatibility: GE DINAMAP  
**Part # 300619 (Qty: 2)**



**5/32" ID Barbed to  
Male Slip Luer, Metal**  
Compatibility: Female Locking Luer  
and Screw Connectors  
**Part # 300667 (Qty: 10)**

## Wash Plugs and O-Rings

GE Healthcare provides wash plugs that block water from entering the BP cuff during cleaning.



**Washing Plugs for  
Submin Connector Cuffs**  
Compatibility: GE,  
GE Marquette, Protocol  
**Part # 330072 (Qty: 50)**



**Washing Plugs for  
Screw Connector Cuffs**  
Compatibility: GE DINAMAP  
**Part # 300877 (Qty: 50)**



**O-Ring Replacement Set for  
Air Hose Instrument Connector**  
Compatibility: (330073) and Air Hoses  
(107365, 107368, 107366)  
**Part # 330077 (Qty: 20)**



**Washing Plugs for Neo-Snap  
Female Recept and DINACLICK  
Male Plug**  
Compatibility: GE  
**Part # 2059295-001 (Qty: 50)**

## About GE Healthcare

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter – great people and technologies taking on tough challenges. From medical imaging, software & IT, patient monitoring and diagnostics to drug discovery, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

GE Healthcare  
8200 West Tower Avenue  
Milwaukee, WI 53223  
U.S.A.

GE Healthcare Finland Oy  
Kuortaneenkatu 2  
00510 Helsinki, Finland  
Europe

GE (China) Co. Ltd - Healthcare  
Room 1105  
Shanghai MAXDO Center  
No. 8 Xing Yi Road  
Hong Qiao Development Zone  
Shanghai, China 200336  
Asia

## For More Information

To learn more about GE Healthcare's comprehensive cuff solutions, visit us at [cuffs.gehealthcare.com](http://cuffs.gehealthcare.com)

© 2013 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

\*Trademarks of General Electric company.

Colin is a registered trademark of Colin Medical Instruments.

Datascope is a registered trademark of Datascope Corp.

HP is a registered trademark of Agilent Technologies.

Nihon Kohden is a registered trademark of Nihon Kohden Corporation.

Philips is a registered trademark of Royal Philips Electronics.

Protocol is a registered trademark of Welch Allyn, Inc.

Siemens is a registered trademark of Siemens Medical Solutions.

Spacelabs is a registered trademark of Spacelabs Medical Systems.

GE Healthcare, a division of General Electric Company.

GE Medical Systems *Information Technologies*, Inc., doing business as GE Healthcare.



imagination at work